

Proximal Lateral Tibia Locking Plate



Bonetech Medisys
implanting precision



Bonetech Medisys
implanting precision

Proximal Lateral Tibia Locking Plate

Two 2.0 mm holes for K-wire and sutures

Anatomically precontoured low profile plates

Combi holes combine a DCU hole with a threaded locking hole

Tapered tip

Indication :

- Proximal shaft fractures
- Metaphyseal fractures
- Intra-articular fractures
- Periprosthetic fractures

HOLE	CODE (Left)	CODE (Right)
5	BT353.T05L	BT353.T05R
7	BT353.T07L	BT353.T07R
9	BT353.T09L	BT353.T09R
11	BT353.T11L	BT353.T11R
13	BT353.T13L	BT353.T13R
15	BT353.T15L	BT353.T15R
T - Titanium, S - Steel		