

Auxein

2.4MM WISE-LOCK MODULAR HAND & FOOT SYSTEM



NEXT-LEVEL CARE

ABOUT US

Auxein Medical is an integrated, research based, orthopaedic implants & instruments manufacturing company, producing a wide range of quality, affordable generic implants, trusted by healthcare professionals and patients across geographies. It is the Company's constant endeavour to provide a wide basket of generic and our innovator products that exceed the highest expectations of customers in term of quality and safety. The company has world-class manufacturing unit established in India and serves customers in over 75 countries worldwide.



IMPACTING HEALTHCARE GLOBALLY

Our Reach & Influence

500+

Employees Worldwide

75+

Operations in Countries

150+

Distributors Worldwide

20_{M+}
Patients Cured

200+

FDA Approved Products

250+

Machines Worldwide

3000+

Products Orthopedic Solutions

Auxein is dedicated to advancing medical technology through innovative orthopedic solutions, including **Trauma, Arthroscopy, and Reconstructive Procedures**, with a commitment to improving patient outcomes worldwide. Our approach integrates advanced research, precision engineering, and high-quality manufacturing to develop reliable implants and surgical instruments that address the evolving needs of modern orthopedic practice.


OUR ACHIEVEMENTS




OVER THE YEARS

Disclaimer: Not all products in this catalogue are covered by European CE, US FDA, or other certifications. Please contact our Sales Team for product-specific certification details.





At Auxein, it is our goal to provide the best in class quality and value to every individual. Consistency is a virtue we pursue with a passion. From the procurement of raw materials to the manufacturing of our range of 1200+ orthopaedic implant products, we employee world-class systems and technology at par with the highest standards of globally recognized evaluation systems. The 2.4mm WISE-LOCK MODULAR HAND & FOOT SYSTEM is part of the orthopedic trauma portfolio, designed for fixation of fractures, osteotomies, and arthrodesis procedures in the small bones of the hand and foot. Developed with clinical input, the system provides stable fixation, controlled load distribution, and supports restoration of anatomical alignment.



WISELOCK

2.4mm Wise-Lock Modular
Hand System





2.4mm Wise-Lock Plate, Straight

STERILE



Holes	Sterile	Non-Sterile
4	TI-505.304-S	TI-505.304
5	TI-505.305-S	TI-505.305
6	TI-505.306-S	TI-505.306
7	TI-505.307-S	TI-505.307
8	TI-505.308-S	TI-505.308
10	TI-505.310-S	TI-505.310



2.4mm Adaption Plate, Straight

STERILE



Holes	Sterile	Non-Sterile
12	TI-821.312-S	TI-821.312




2.4mm Wise-Lock Adaption Plate, Straight

STERILE



Holes	Sterile	Non-Sterile
12	TI-555.312-S	TI-555.312



2.0mm Wise-Lock Condylar Plate

STERILE



Holes	Sterile	Non-Sterile
Shaft 4, Head 2	TI-559.304-S	TI-559.304
Shaft 5, Head 2	TI-559.305-S	TI-559.305
Shaft 6, Head 2	TI-559.306-S	TI-559.306
Shaft 7, Head 2	TI-559.307-S	TI-559.307



Note: Use corresponding screws as per the specific plate profile and system compatibility.



2.4mm Wise-Lock T-Adaption Plate

STERILE



Holes	Sterile	Non-Sterile
Shaft 7, Head 2	TI-555.307T-S	TI-555.307T



2.4mm Wise-Lock T-Plate, Titanium

STERILE



Holes	Sterile	Non-Sterile
Shaft 4, Head 3	TI-556.304-S	TI-556.304
Shaft 5, Head 3	TI-556.305-S	TI-556.305
Shaft 6, Head 3	TI-556.306-S	TI-556.306
Shaft 7, Head 3	TI-556.307-S	TI-556.307



2.4mm Wise-Lock Y-Adaption Plate


STERILE



Holes	Sterile	Non-Sterile
Shaft 4, Head 3	TI-555.304Y-S	TI-555.304Y
Shaft 5, Head 3	TI-555.305Y-S	TI-555.305Y
Shaft 6, Head 3	TI-555.306Y-S	TI-555.306Y
Shaft 7, Head 3	TI-555.307Y-S	TI-555.307Y



Note: Use corresponding screws as per the specific plate profile and system compatibility.

The background of the page features a light gray, wavy, topographic-like pattern that flows from the top and bottom edges towards the center, framing the central text.

SCREWS



2.4mm Cortical Screw, Self-Tapping, (Star Head)



Length	Code
6mm	TI-1144.06
7mm	TI-1144.07
8mm	TI-1144.08
9mm	TI-1144.09
10mm	TI-1144.10
11mm	TI-1144.11
12mm	TI-1144.12
13mm	TI-1144.13
14mm	TI-1144.14
16mm	TI-1144.16
18mm	TI-1144.18
20mm	TI-1144.20
22mm	TI-1144.22
24mm	TI-1144.24
26mm	TI-1144.26
28mm	TI-1144.28
30mm	TI-1144.30
32mm	TI-1144.32
34mm	TI-1144.34
36mm	TI-1144.36
38mm	TI-1144.38
40mm	TI-1144.40

STERILE

Length	Code
6mm	TI-1144.06-S
7mm	TI-1144.07-S
8mm	TI-1144.08-S
9mm	TI-1144.09-S
10mm	TI-1144.10-S
11mm	TI-1144.11-S
12mm	TI-1144.12-S
13mm	TI-1144.13-S
14mm	TI-1144.14-S
16mm	TI-1144.16-S
18mm	TI-1144.18-S
20mm	TI-1144.20-S
22mm	TI-1144.22-S
24mm	TI-1144.24-S
26mm	TI-1144.26-S
28mm	TI-1144.28-S
30mm	TI-1144.30-S
32mm	TI-1144.32-S
34mm	TI-1144.34-S
36mm	TI-1144.36-S
38mm	TI-1144.38-S
40mm	TI-1144.40-S





2.4mm Wise-Lock Screw, Self-Tapping, (Star Head)




Length	Code
6mm	TI-1145.06
7mm	TI-1145.07
8mm	TI-1145.08
9mm	TI-1145.09
10mm	TI-1145.10
11mm	TI-1145.11
12mm	TI-1145.12
13mm	TI-1145.13
14mm	TI-1145.14
16mm	TI-1145.16
18mm	TI-1145.18
20mm	TI-1145.20
22mm	TI-1145.22
24mm	TI-1145.24
26mm	TI-1145.26
28mm	TI-1145.28
30mm	TI-1145.30
32mm	TI-1145.32
34mm	TI-1145.34
36mm	TI-1145.36
38mm	TI-1145.38
40mm	TI-1145.40
42mm	TI-1145.42
44mm	TI-1145.44
46mm	TI-1145.46
48mm	TI-1145.48
50mm	TI-1145.50

STERILE

Length	Code
6mm	TI-1145.06-S
7mm	TI-1145.07-S
8mm	TI-1145.08-S
9mm	TI-1145.09-S
10mm	TI-1145.10-S
11mm	TI-1145.11-S
12mm	TI-1145.12-S
13mm	TI-1145.13-S
14mm	TI-1145.14-S
16mm	TI-1145.16-S
18mm	TI-1145.18-S
20mm	TI-1145.20-S
22mm	TI-1145.22-S
24mm	TI-1145.24-S
26mm	TI-1145.26-S
28mm	TI-1145.28-S
30mm	TI-1145.30-S
32mm	TI-1145.32-S
34mm	TI-1145.34-S
36mm	TI-1145.36-S
38mm	TI-1145.38-S
40mm	TI-1145.40-S
42mm	TI-1145.42-S
44mm	TI-1145.44-S
46mm	TI-1145.46-S
48mm	TI-1145.48-S
50mm	TI-1145.50-S



The background of the page features a light gray, wavy, topographic-like pattern that flows across the top and bottom edges, framing the central text.

INSTRUMENTS




1137-01

Self-Centering Double Drill Guide Ø2.4/1.8mm



1137-02

Wise-Lock Drill Sleeve 2.4mm For Drill Bit Ø1.8mm



1136-06

Countersink For 2.0/2.4/2.7mm Screws



1137-03

Bending Pin With Thread For 2.4mm/2.7mm Wise-Lock Plates



1137-04

Drill Bit Quick Coupling Ø1.8mm x Length 110/82mm, 2- Flute



1137-05

Drill Bit Quick Coupling Ø2.4mm x Length 100/73mm, 2-Flute





1137-07

Screwdriver Shaft, T8 For 2.4mm/2.7mm Screws



1137-08

Holding Sleeve For 2.4mm Cortical Screws




1137-09

Holding Sleeve For 2.4/2.7mm Wise-Lock Screws



1138-10

Container with Tray and Lid For 2.4mm Implant & Instrument Set





2.4MM WISE-LOCK MODULAR HAND & FOOT SYSTEM

The content of this catalogue is intended for healthcare professionals only. All product information, illustrations, descriptions, and specifications provided here are for general informational purposes and are subject to change without notice. Auxein Medical Pvt. Ltd. reserves the right to improve or discontinue products at any time without prior notification.

This document does not constitute medical advice. Surgeons and medical practitioners must exercise their own professional judgment when using any medical devices or instruments, including those provided by Auxein.

INTENDED USE

The implants and instruments in this catalogue are intended for use by qualified orthopedic surgeons trained in fracture fixation procedures. Proper implant selection and surgical technique are essential for safe and effective use. Improper selection, misuse, or insufficient surgical experience may result in complications.

REGULATORY COMPLIANCE CERTIFIED

CE Approved (EUMDR-2017/745)

USFDA 510(k)

ISO 13485:2016


CDSCO

MDSAP

DISCLAIMER: Not all products in this catalogue are covered by European CE, US FDA, or other certifications. Please contact our Sales Team for product-specific certification details.

INTELLECTUAL PROPERTY

All trademarks, trade names, brand names, product names, and logos mentioned in this catalogue are the property of Auxein Medical Pvt. Ltd. Unauthorized use or reproduction is strictly prohibited.



SCAN TO KNOW MORE

INDIA
AUXEIN MEDICAL PVT.LTD.

PLOT NO. 168-169-170, PHASE-4,
KUNDLI INDUSTRIAL AREA,
HSIIDC, SECTOR-57,
SONEPAT-131028, HARYANA
TEL: +91 99106 43638
FAX: +91 86077 70197
EMAIL: INFO@AUXEIN.COM

MÉXICO
AUXEIN MÉXICO, S.A DE CV

TEPIC 39 INT 801, COLONIA
ROMA SUR,
ALCADÍA CUAUHTÉMOC, CDMX
MÉXICO, C.P 06760
TEL: +521 55 7261 0318
EMAIL: INFO@AUXEIN.MX

USA
AUXEIN INC.

1500 NW 89TH COURT,
SUITE 107-108
DORAL, FLORIDA 33172
TEL. +1 305 395 6062
E FAX: +1 305 395 62 62
EMAIL: USOFFICE@AUXEIN.COM

DUBAI UAE
AUXEIN MEDICAL

TRADING FZCO OUD METHA
OFFICES BUILDING B BLOCK,
1ST FLOOR, OFFICE NO.26,
NEAR HEALTHCARE CITY,
PO BOX-25572
TEL. +971 438 71333
EMAIL: INFO@AUXEIN.COM

WWW.AUXEIN.COM